COMPRESSOR NEBULIZER

MODEL TURBO BASE
USER'S MANUAL

- EN COMPRESSOR NEBULIZER**
Model TURBO BASE
User's manual (3-7)
- DE DRUCKZERSTÄUBER**
Modell TURBO BASE
Benutzerhandbuch (7-12)
- FR NÉBULISEUR À COMPRESSEUR**
Modèle TURBO BASE
Manuel d'utilisation (12-17)
- IT NEBULIZZATORE A COMPRESSORE**
Modello TURBO BASE
Manuale illustrativo (17-22)
- RO NEBULIZATOR CU COMPRESOR**
Model TURBO BASE
Manual de utilizare (22-27)
- MD**
- BG КОМПРЕСОРЕН ИНХАЛАТОР**
Модел TURBO BASE
Инструкция за употреба (27-32)
- PL NEBULIZATOR SPRĘŻARKOWY**
Model TURBO BASE
Instrukcja obsługi (32-37)

EN DEAR DR. FREI COMPRESSOR NEBULIZER OWNER

Thank you for choosing Dr. Frei® Compressor Nebulizer model TURBO BASE. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Swiss trademark Dr. Frei®. Your nebulizer is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy. Aerosol easily and quickly comes to both upper and lower respiratory tracts. Before starting to use the compressor nebulizer, please, study the user's manual carefully. Following the rules stated in the user's manual will provide for the faultless operation of the compressor nebulizer for a long period of time. For all questions concerning the unit, please, contact service centre of distributor of trademark Dr. Frei® in your country.

TABLE OF CONTENTS

| | |
|---|---|
| INTRODUCTION..... | 3 |
| What is a nebulizer and how does it work?..... | 3 |
| What is nebulizing therapy used for?..... | 3 |
| What are the advantages of nebulizing therapy?..... | 3 |
| BEFORE STARTING TO USE THE UNIT..... | 3 |
| Important Safety Information..... | 3 |
| Unit Accessories..... | 4 |
| CARRYING OUT NEBULIZING THERAPY..... | 4 |
| Operating the Unit..... | 4 |
| Changing the Air Filter..... | 5 |
| UNIT MAINTENANCE..... | 5 |
| Cleaning the Unit..... | 5 |
| Troubleshooting..... | 5 |
| TECHNICAL SPECIFICATIONS..... | 5 |
| SYMBOL DESCRIPTION..... | 6 |
| WARRANTY..... | 7 |
| ELECTROMAGNETIC COMPATIBILITY INFORMATION..... | 7 |

INTRODUCTION

What is a nebulizer and how does it work?

The compressor nebulizer system is a device that turns a liquid drug into a fine aerosol under the influence of the compressed air. Aerosol contains large and small particles. When breathing the gravity makes the larger particles accumulate on the walls of the upper respiratory passages while the smaller particles have enough time to reach the lower bronchial branches. Therefore the drug comes to all parts of the patient's bronchial tree, including the smallest bronchi and alveoli, and gets into the blood vessels on the walls of the alveoli in big therapeutic doses with tiny or none side effects.

What is nebulizing therapy used for?

Nebulizer systems are used by the patients with:

- acute respiratory diseases,
- obstructive pulmonary disease,
- chronic inflammatory processes in respiratory passages such as chronic bronchitis, bronchial asthma (incl. chronic or unstable ones, exacerbation of asthma), chronic pharyngitis,
- respiratory passages infections,
- abnormalities in secretion quality and secreting in bronchi,
- different types of coughing.

What are the advantages of nebulizing therapy?

Nebulizing therapy gives a possibility to quickly effect the inflammatory zone by big drug doses which increases the influence of the treatment and reduces the risk of side effects to other organs. Nebulizing therapy is the easiest way of inhalation treatment and it is also an economical means of treatments from the point of view of purposeful use of drug. Aerosol treatment with the help of nebulizing systems can be used for children, elderly and infirm patients.

BEFORE STARTING TO USE THE UNIT

Important Safety Information

- Carefully read the instructions before use. Keep this instruction manual for future references.
- Use the device only as inhalator for therapeutic purposes. Medicine type, mode and time of application should be prescribed by your doctor. Installation operations have to be done according to the instruction in this handbook and only for the uses declared from the constructor.
- This device is unsuitable for use with inflammable anaesthetic mixtures containing air, oxygen or nitrogen protoxide.
- Always disconnect the device from the main socket after use and before cleaning operations.

- Do not touch the device if it accidentally falls into the water. Unplug it immediately and do not use afterward. Contact one of authorized service centers.
- Do not touch the device with wet hands.
- During use, make sure that the unit is placed on a flat and stable surface.
- Do not use the device when patient is sleeping.
- Refer all servicing to the service centre. Do not open the device. If the device does not work properly turn it off and consult the user's manual.
- Do not leave any packaging materials (plastic bags, cardboard boxes, and so on) within the reach of children to prevent dangers.
- Keep accessories out of reach of children. Children and people with reduced physical or mental capacities need to use the device only under close surveillance of an adult who has read this instruction manual. Keep the medicine cup out of the reach of children under 36 months because it contains small parts that could be ingested.
- Always disconnect the device from the main socket when it is not used.
- It is recommended to use one personal medication set of accessories.
- Do not leave pharmaceutical residual into the medicine cup and/or into the accessories at the end of the therapy.
- Cleaning and / or maintenance must be done only after turning off the device and unplugged power cord from the outlet.
- Do not pull cable or appliance in order to remove the plug.
- This is a medical device for home use and should be used on medical prescription. Must be run as indicated on this user's instructions. It's important that patients read and understand the information for use and maintenance of the unit. Contact your dealer or service centre for any questions.
- If the equipment's plug does not fit the socket, contact a qualified technician to substitute the plug. Do not use multiple sockets and/or extensions.
- Do not obstruct the air intake.
- The installation must be done according to the instructions and use declared by the manufacturer. Improper installation and / or improper use may cause damage to persons, animals or things. The manufacturer can not be held responsible.
- Do not use the device while bathing and showering.
- Do not expose the device to the action of dust, heat (that could deform the housings), sunlight or other atmospheric agents that could compromise the safety of the device itself.
- Do not use the device nearby sources of potential interference (example: heat sources like fireplaces, electric or gas heaters; steam sources like boilers or kettles, electromagnetic radiation sources like microwave ovens, Wi-Fi devices).
- Keep away from pets that could contaminate the surface in contact with the patient or in general hinder the right performance of the therapy.
- The appliance generates compressed air. Secondary consequences could come for the intolerance to the medicine used.


Electric devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

KEEP THIS INSTRUCTION MANUAL DURING THE WHOLE SERVICE LIFE OF THE UNIT

Unit Accessories

- 1 Mouthpiece
- 2 Adult mask
- 3 Pediatric mask
- 4 Medicine cup - upper part
- 5 Medicine cup - lower part
- 6 Medicine cup - air entrance
- 7 Nozzle
- 8 Air tube
- 9 Spare air filters (5 pcs)


CARRYING OUT NEBULIZING THERAPY

The appliance is a medical device of class Ia, which converts medicine from the fluid form to vapour in order to adBASEster it directly into the respiratory tracts. The appliance nebulizes medicines both in suspension and in solution. The therapeutic effect could be reduced with oily and/or pasty medicines. The appliance is NOT a nebulizer for essential oil. Don't preheat liquid and/or alcohol solutions. Medicine cup and all accessories have NOT been sterilized. Read «Unit maintenance» paragraph before use.

Operating the Unit

1. Put the device on a flat, stable surface.
2. Ensure that the rated voltage corresponds to the mains voltage supply.
3. Connect the device into the main socket.
4. Open the medicine cup unscrewing the two parts (figure A).
5. Insert the nozzle in the LOWER PART OF THE medicine cup (figure B).
6. Fill the medicine cup with pharmaceutical prescribed by your doctor (figure C). Note: The notches on the lower part of medicine cup are only indicative.
7. Close the medicine cup screwing two parts back (figure D).
8. Connect one end of the air tube (9) to the medication cup and the other end to the connector on the device (Figure E).
9. Install the needed accessory on the medicine cup: adult or pediatric mask, or mouthpiece (figure F).
10. Turn on the device moving the switch to ON position.
11. After the application, turn off the device and unplug. Before restarting the appliance make sure it has been cooled to room temperature.


RECOMMENDATIONS FOR CARRYING OUT NEBULIZING THERAPY:

- Sit down in a comfortable and upright position.
- If a mask (adult or pediatric) is used, be sure that it securely adheres to the face. Mask has holes for the air flow during the expiration.
- Breathe easy and have a short pause at the end of each breath.
- Therapeutic application ends when the medicine came out in a discontinuous mode (spluttering).
- When using the mouthpiece, do not tilt the medicine cup over 45°.

Changing the Air Filter

The filter has to be substituted after 30 hours of use or when it becomes grey. Spair air filter and its container are on the figure H. For filter replacing, carefully raise the lid using a flat screwdriver as shown in the picture I. Remove the filter and fit in a new filter (9). Fit the filter carrier into its housing making sure it is attached securely.

UNIT MAINTENANCE

Cleaning the Unit

After each therapy remove the device from the socket and clean the device and the accessories. Attention: use only the accessories supplied with the product, other accessories can lead to ineffective therapy.

CLEANING OF THE MOUTHPIECE

Use drinking water to carefully clean these accessories. They can be sterilized by boiling for 5-10 minutes into the water or by immersion in a chemical disinfectant (oxidizing electrolytic chloride), following product instructions and rinsing off with drinking water. Leave to dry naturally in a clean place.

CLEANING OF THE MEDICINE CUP

After use, open the two parts, wash with soap and hot water and rinse it carefully. Remove excess of water with a soft, clean and dry cloth and leave to dry naturally in a clean place. To sterilize, cover the medication container with a cold sterilizing solution, following the instructions on the sterilizer pack. It is recommended to use one personal medication container for each patient, to be use for 6 months or 120 treatments. Do not use boiling water to clean it. Carefully wash masks and hoses with disinfectant solution only.

CONTAMINATION BY MICROBES

In case of pathologies of microbial or viral origin, where infection is possible, the device and its accessories should not be shared with other patients.

Troubleshooting

| PROBLEM | ACTION |
|--|---|
| The device is not turning on. | <ul style="list-style-type: none">• Make sure that the plug has been correctly inserted in the main supply.• Make sure that the switch is in "ON" position. |
| The device is turned on but it's not nebulizing. | <ul style="list-style-type: none">• Make sure to have inserted the nozzle (7) into the medicine cup.• Make sure that the air tube is not bended or crushed (8).• Make sure that the air intake is not obstructed and the air filter is not dirty. Replace the air filter with the new one if necessary (9).• Make sure that the medicine cup has been filled in with medicine. |

NOTE: If the device is not functioning correctly in spite of inspections effected, contact the TM Dr. Freil® service center in your country.

TECHNICAL SPECIFICATIONS

Compressor nebulizer machine. Complies with the requirements of European Standard EN 13544-1. Piston rotary presser, without lubrication. Medical device of class IIA according to Directive 93/42/EEC. Medicine Linear Cup.

| | |
|-----------------|--|
| Model | Turbo Base |
| Adapter input | 100 - 240 V ~ , 50 / 60 Hz, max 0.5 A |
| Adapter output | 12.0 V  1 A |
| Main unit Input | DC 12V, 1A |
| Noise Level | ≤ 55 dB |


| | |
|---|--|
| Max capacity of nebulizer kit | 10 ml |
| Particle size MMAD | 3 μm |
| Nebulizing rate | ~ 0.3 ml/min |
| Operation mode | 20 minutes on, 40 minutes off |
| Operation temperature, humidity and atmospheric pressure | 10°C ~ 40°C, $\geq 85\%$ relative humidity 860 ~ 1060 hPa |
| Transport and storage temperature humidity and atmospheric pressure | -10°C ~ 40°C, $\geq 85\%$ relative humidity 860 ~ 1060 hPa |
| Pollution degrees | Degrees 2 |
| Overvoltage category | Category II |
| Maximum operational altitude, m | 2000 |
| Accessories | Compressor nebulizer machine, medicine cup, nozzle, adult mask, pediatric mask, mouthpiece, air tube, AC adapter |
| Size, mm | 110 (W) x 110 (H) x 62 (D) |

* Subject to technical modification without prior notice.

This device has an expected service life of about 2,000 therapy cycles, basing on an average therapy duration of 8 min.

Electromedical devices require special care. During installation and use with respect to EMC requirements, it therefore required that they are installed and/or used according to the manufacturer's specification. There is a potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electro-medical device.

SYMBOL DESCRIPTION

| SYMBOL | REFERENT | SYMBOL | REFERENT |
|---|---|---|---|
|  | European representative |  | CE mark |
|  | Double insulation |  | Do not use the appliance when taking a bath or shower |
|  | BF applied part | ON | Device is ON |
|  | Obligation to read instruction manual | OFF | Device is OFF |
|  | General warning symbol |  | Consult instruction for use |
| IP 21 | Protection against solid and liquid bodies penetration (protected against solid bodies over 12 mm; protected against access with a finger; protected against vertically falling water drops). |  | Keep dry |
|  | This product (including its related accessories) must be delivered at the end of its service life at a collection point for the recycling by electric and electronic devices, by the user or by His dealer. |  | Temperature limitation |
|  | Manufacturer |  | Manufacturing date |

| | | | |
|---|----------------|--|---------------------|
|  SN | Serial number |  | Alternating current |
|  | Direct current | | |

WARRANTY

This device is guaranteed for 2 years from the date of purchase. The warranty does not cover consumable materials, and those parts subject to normal wear and tear. The warranty is only valid upon presentation of the warranty card with purchase date and stamp of trade organization and receipt.

Should the unit be damaged due to user misuse or negligence, the manufacturer or authorized dealers shall not be responsible for any loss or damage to the unit.

All assembly, extension, adjustment or repair of this unit must be performed only by the authorized service center of **Dr. Frei®** distributor in your country.

For repairs or spare parts purchase, please, contact the authorized service center of **Dr. Frei®** distributor in your country.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

This unit conforms to EMC Standard IEC60601-1-2.

Guidance and manufacture's declaration – electromagnetic emission

The Turbo Base is intended for use in the electromagnetic environment specified below. The user of the Turbo Base should assure that it is used in such environment.

| Emission test | Compliance | Electromagnetic environment – guidance |
|---|------------|---|
| RF emissions CISPR 11 | Group 1 | The Turbo Base uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emission CISPR 11 | Class B | The Turbo Base is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies | |

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■ HONSUN (NANTONG) Co., Ltd. No.8 Tongxing Road, Economic&Technological, Development Area,
☐ Shanghai International Holding Corp. GmbH (Europe), Eiffelstraße 80, 20537 Hamburg, Germany.